



Connecting innovation.
Delivering better outcomes.



Shield Style Goggles

Fits over prescription eyewear. Wrap around shield with integrated side protection.


ADDRESS 36/2 Chaplin Drive, Lane Cove NSW 2066, Australia
PHONE +61 2 8313 5478
EMAIL info@innoviagroup.com.au
WEB innoviagroup.com.au
ABN 56 688 008 380

LASER SAFETY EYEWEAR

THULIUM AND HOLMIUM LASER PROTECTION GOGGLES

THIS shield style eyewear is for 1900 to 3000nm laser use. TR90 frame with polycarbonate filter, 55% visible light transmittance, and OD5+ across the 2000 to 2300nm band used by thulium and holmium systems.

PROTECTION RANGE	1900 to 3000nm
FRAME MATERIAL	TR90 frame, PC filter
VISIBLE LIGHT TRANSMITTANCE	55% VLT
OPTICAL DENSITY	1900 to 3000nm OD4+ · 2000 to 2300nm OD5+
LB RATING	1900 to 3000nm DI LB3 LP S CE
APPLICATION	Thulium laser, holmium laser



Innovia Group Pty Ltd · Medical device distribution, Australia and New Zealand
Specifications from the manufacturer filter datasheet. Full test data available on request.

innoviagroup.com.au